

K-ASSAY®

hsCRP Control

Lot 1200, Exp. 2022-04-27

Cat. No. K180C-2M**INTENDED USE**

The **K-ASSAY®** hsCRP Control is intended for use as an assayed quality control material for monitoring the performance of C-Reactive Protein (CRP) immunoturbidimetric assays. *FOR IN VITRO* DIAGNOSTIC USE.

SUMMARY

The use of quality control materials to objectively monitor the precision of procedures in use in the clinical laboratory has been well established. The **K-ASSAY®** hsCRP Control is provided at two levels to assist in the monitoring of analytical systems within the clinical range.

SET COMPOSITION

K180C-2M
Level 1, 2 Human serum (liquid) 1 x 3 mL each level

hsCRP Control Levels 1 and 2 contain pooled human serum with assigned values for CRP.

WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY. Rx only.

Not to be used internally in humans or animals. Normal precautions exercised in handling laboratory reagents should be followed.

Controls contain pooled human serum from CRP positive human serum. The serum has been tested and found non-reactive for the presence of HBsAg, HCV Ab, and HIV antibody by an FDA approved method. However, it is not possible to guarantee that any human source material is free of these infectious agents.

Therefore, all products containing human serum should be handled in accordance with good laboratory practices using appropriate precautions as described in the CDC/NIH Manual, "Biosafety in Microbiological and Biomedical Laboratories," 2nd ed., 1988, HHS Publication No. (CDC) 88-8395.

Do not mix or use controls from one test set with those from a different lot number.

Do not use controls past their expiration date stated on each control container label.

Do not pipette by mouth. Avoid ingestion and contact with skin.

Controls in this set contain < 0.1% w/v sodium azide as a preservative. Sodium azide may form explosive compounds in metal drain lines. When disposing of controls through plumbing fixtures, flush with copious amounts of water. For

further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control, Atlanta, GA.

CONTROL PREPARATION

1. Gently invert the bottle before use.
2. The control is ready to use.
3. Replace the cap immediately after use to prevent evaporation or contamination.

STORAGE AND HANDLING

Store controls at 2-8°C. Return all controls to 2-8°C promptly after use. Controls can be used up to the expiration date shown on the vial label and package insert.

CONTROL STABILITY

Discard controls if they become contaminated. Evidence of cloudiness or particulate material in solution is cause to discard.

PROCEDURE**Materials Supplied**

Level 1	Human serum	1 x 3 mL
Level 2	Human serum	1 x 3 mL

Materials Required But Not Supplied

K-ASSAY® hsCRP Assay

K-ASSAY® hsCRP Calibrator

or

K-ASSAY® CRP (3) Assay

K-ASSAY® CRP (3) Calibrator

Two-reagent clinical chemistry analyzer

Assay Procedure:

NOTE: Allow all reagents and specimens to come to room temperature. Mix all reagents gently before using.

K-ASSAY® hsCRP Controls are assayed using the same procedure as the patient test samples run in the test procedure. See package insert from the **K-ASSAY®** hsCRP or **K-ASSAY®** CRP (3) immunoturbidimetric assay.

LIMITATIONS OF PROCEDURE

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, and good laboratory technique. Erroneous results can occur from improper storage, reconstitution inaccuracies, and technical errors associated with assay procedures.

This product is intended for use as an assayed control for quantitative assays of listed constituents in human serum or plasma. This product is not intended for use as a calibrator.

EXPECTED VALUES

Actual values recovered depend on the instrument and reagent used. The assignment of mean values was derived from analysis of vials representative of the entire lot.

Values listed are specific for each lot only. Verify that the lot numbers on the vials of hsCRP Control correspond to the lot numbers listed for the Assay Data.

The Expected Range and the Mean are provided to assist the laboratory until it has established its own mean and standard deviation. It is considered good laboratory practice for each laboratory to establish its own mean and standard deviation for its test methods. The indicated Mean and Expected Range of the Mean should serve as a guide in assessing the performance of each test method.

Recovered values may be method dependent. The variations that can occur over time and between laboratories may be attributed to differences in laboratory technique, instrumentation, reagent lot, method modifications, and other systematic errors including random errors.











ASSAY DATA

Lot : 1200

L1		L2		UNITS
MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE	
0.211	0.169 - 0.253	0.807	0.686 - 0.928	mg/dL
2.11	1.69 - 2.53	8.07	6.86 - 9.28	mg/L
211	169 - 253	807	686 - 928	µg/dL

The expected values for the **K-ASSAY®** hsCRP Control are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert for each lot for the appropriate control values.

LABELING SYMBOLS

	Lot Number
	Control
	Expiration or "Use By" Date
	Catalog Number
	For <i>In Vitro</i> Diagnostic Use
	2-8°C Temperature Limitation. Store between 2 and 8 degrees C
	Potential Human Biohazard
	Manufacturer
	Consult Package Insert for Instructions for Use
	Authorized Representative in the European Community

EU AUTHORIZED REPRESENTATIVE





Advena Ltd.

Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta

ORDERING / PRICING / TECHNICAL INFORMATION



KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive
Seattle, WA 98168 USA
TEL: (206) 575-8068 / (800) 526-4925
FAX: (206) 575-8094

Printed February 2021