



**KAMIYA BIOMEDICAL COMPANY**

# Nitric Oxide Fluorometric Assay

**For the quantitative determination of total nitric oxide (NO) in culture media,  
plasma and tissue homogenates**

**Cat. No. KT-152**

**For Research Use Only. Not For Use in Diagnostic Procedures.**

## PRODUCT INFORMATION

### Nitric Oxide Fluorometric Assay

Cat. No. KT-152

#### PRODUCT

The **K-ASSAY®** Nitric Oxide Fluorometric Assay is for the quantitative determination of total nitric oxide (NO) in culture media, plasma and tissue homogenates.

#### PRINCIPLE

Nitric oxide (NO) plays an important role in neurotransmission, vascular regulation, immune response and apoptosis. Since NO is rapidly converted to nitrite ( $\text{NO}_2^-$ ) and nitrate ( $\text{NO}_3^-$ ), the total concentration of nitrite and nitrate is used as a quantitative measure of NO production. The **K-ASSAY®** Nitric Oxide Fluorometric Assay provides an accurate and convenient measurement of total nitrate/nitrite concentration in a simple two-step process. In the first step, nitrate is converted to nitrite by nitrate reductase. In the second step, nitrite reacts with the fluorescent probe DAN (2, 3-diaminonaphthalene). NaOH enhances the fluorescent yield. The fluorescent intensity is proportional to the total nitric oxide production. The kit has been tested with culture media, plasma and tissue homogenates.

#### COMPONENT

• Assay Buffer	1 bottle	30 mL
• Enzyme Cofactors (blue cap)	1 vial	Lyophilized
• Enhancer (purple cap)	1 vial	Lyophilized
• Nitrate Reductase (green cap)	1 vial	Lyophilized
• Nitrate Calibrator (yellow cap)	1 vial	Lyophilized
• Nitrite Calibrator (orange cap)	1 vial	Lyophilized
• DAN Reagent (red cap)	1 vial	1 mL
• Sodium Hydroxide (clear cap)	1 vial	1 mL
• Microtiter Plate	2 plates	
• Plate Cover	2 covers	

#### PROTOCOLS

##### Preparation of Reagents

1. Assay Buffer: The assay buffer is ready to use as supplied. Store at 4°C.
2. Enzyme Cofactors: Reconstitute with 1.2 mL of Assay Buffer. Aliquot desired amount and store at -20°C. Keep on ice during use. Freeze/thaw should be limited to 1 time.
3. Enhancer: Reconstitute with 1.2 mL of Assay Buffer. Keep on ice during use. Store at -20°C.
4. Nitrate Reductase: Reconstitute with 1.2 mL of Assay Buffer. Aliquot desired amount and store at -20°C. Keep on ice during use. Freeze/thaw should be limited to 1 time.
5. Nitrate/Nitrite Calibrators: Reconstitute with 1.0 mL of Assay Buffer. Vortex to generate 10 mM calibrator each. Store at 4°C when not in use (do not freeze!). The reconstituted calibrators are stable for 4 months when stored at 4°C.
6. Fluorometric DAN Reagent and NaOH: Ready to use. Store at 4°C.

##### Measurement of Nitrate + Nitrite:

1. Prepare Calibrators:
  - Add 5  $\mu\text{L}$  of the reconstituted 10 mM nitrate/nitrite calibrators to 995  $\mu\text{L}$  assay buffer.
  - Vortex to generate 50  $\mu\text{M}$  working calibrator solution.
  - Add 0, 4, 8, 12, 16, 20  $\mu\text{L}$  of the working calibrator to 6 consecutive wells to generate 0, 200, 400, 600, 800, 1,000 pmoL/well calibrator.
  - Bring the volume to 75  $\mu\text{L}$  with Assay Buffer.

\*Note: DAN Reagent reacts with nitrite, not nitrate. For routine total nitrite/nitrate assay, you may prepare a nitrate calibration curve only. However, if you need to measure nitrite and nitrate concentrations separately, you may prepare a nitrite calibration curve in the absence of Nitrate Reductase in the calibration curve and assay samples. Nitrate = Total – Nitrite.

2. Prepare Samples:

- Samples containing high protein concentration may need to be filtered through a 10Kd MW cut-off filter prior to assay.
- Add 0-75  $\mu\text{L}$  of sample to the wells and adjust the volume to 75  $\mu\text{L}$  with Assay Buffer.

\*Note: Typical urine nitrite and nitrate levels are in the 0.2-2 mM and 1-20  $\mu\text{M}$  range respectively. Typical normal serum levels are 0-20  $\mu\text{M}$  and 0-2  $\mu\text{M}$  respectively with various disease states elevating these levels significantly. Plasma samples or tissue homogenates should be assayed with no more than 10  $\mu\text{L}$  of undiluted sample.

3. Add 5  $\mu\text{L}$  of the Enzyme Cofactor mixture to all wells.
4. Add 5  $\mu\text{L}$  of the Nitrate Reductase to nitrate assay wells (unknowns and calibrators), add 5  $\mu\text{L}$  of buffer in place of Nitrate Reductase (unknowns and calibrators) when you determine nitrite separately.
5. Cover the plate with the plate cover and incubate at room temperature for 1-4 hrs. 1 hour = ~90% conversion of Nitrate to Nitrite, 2 hours = ~95% conversion, 4 hours = ~99% conversion.
6. Add 5  $\mu\text{L}$  of Enhancer to each well. Incubate 30 minutes to quench interfering compounds.
7. Add 5  $\mu\text{L}$  of DAN Reagent to each well. Incubate for 10 minutes at room temperature.
8. Add 5  $\mu\text{L}$  of NaOH to each well. Incubate for 10 minutes at room temperature.
9. Read the plate in a fluorometer using Ex. = 360 nm and Em. = 450 nm.

## CALCULATIONS

1. Plot calibration curve: Plot fluorescence vs. picomoles nitrate.
2. Determine sample nitrate and nitrite concentrations:

$$C = [\text{Nitrate} + \text{Nitrite}] (\mu\text{M}) = \left( \frac{\text{fluorescence} - y \text{ intercept}}{\text{slope}} \right) \left( \frac{1}{\text{sample volume}(\mu\text{L})} \right) \times \text{dilution} *$$

OR

$C = S_a/S_v$ , where  $S_a$  is the amount of samples as read from calibration curve (in pmole), and  $S_v$  is the volume of sample added to the well (in  $\mu\text{L}$ ), multiplied by the dilution\* factor.

\*Dilution is the sample dilution done prior to addition of the sample to the well.

## STORAGE

Store unopened kit at  $-20^\circ\text{C}$ .

### FOR RESEARCH USE ONLY

## KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle WA 98168  
 Tel: (206) 575-8068 Fax: (206) 575-8094  
 Email: LifeScience@k-assay.com  
 www.k-assay.com