

**KAMIYA BIOMEDICAL COMPANY**

# Rat Osteoclast Culture Kit

**For the culture of Osteoclasts from precursor cells.**

**Cat. No.: KT-361, KT-703**

**For Research Use Only.**

**PRODUCT INFORMATION**

**Rat Osteoclast Culture Kit**  
**Cat. No. KT-361, KT-703**

**PRINCIPLE**

In aging societies, osteoporosis and other age-related bone metabolism disorders are a rapidly increasing problem. The amount of bone in an organism is controlled by a balance of osteoblasts (bone-forming cell) and osteoclasts (bone-destroying cell) activities. A method to induce osteoclasts formation from bone marrow cells using M-CSF (macrophage-colony stimulating factor) and RANKL (receptor activator of NF-κB ligand) has been established in recent years. This kit includes cryopreserved primary precursor osteoclasts from rat bone marrow and Culture Medium containing M-CSF and RANKL.

**COMPONENTS**

**Rat Osteoclast Culture Kit**

Components	KT-361	KT-703
Rat Osteoclast Precursor Cells, frozen	2 vials with 2 x 10 <sup>6</sup> cells	4 vials with 2 x 10 <sup>6</sup> cells
Washing Medium, alpha-MEM	50 mL	100 mL
Culture Medium, alpha-MEM with M-CSF (50 ng/mL) and RANK Ligand (15 ng/mL)	25 mL	50 mL

**Storage**

Components	Storage Conditions	Shelf-Life
Rat Osteoclast Precursor Cells	Liquid Nitrogen (preferred)	See expiration date on box
Washing Medium	-20°C Freezer	6 months
	-80°C Freezer	1 year
Culture Medium	-20°C Freezer	6 months
	-80°C Freezer	1 year

**Materials required but not provided**

- Pipettes
- 96-well, flat bottom culture plate
- Tubes
- Refrigerated centrifuge
- Water bath

**PRECAUTIONS**

1. Read the instructions carefully before beginning the culture.
2. This kit is for research use only, not for human or diagnostic use.

**PROTOCOL**

1. Thaw a vial of primary precursor osteoclasts in a 37°C water bath.
2. After thawing, transfer the cells to a 15 mL centrifuge tube, add 10 mL of Wash Medium and mix briefly. Centrifuge 1,000 rpm for 5 minutes at 4°C.

3. Remove supernatant and add 10 mL of Wash Medium and mix briefly. Centrifuge 1,000 rpm for 5 minutes at 4°C.
4. Remove supernatant and resuspend the cells in 2.5 - 5 mL of Culture Medium. To study factors that effect osteoclasts formation, add the factors to the Culture Medium.
5. Transfer 100  $\mu$ L of cell suspension into each well of a 96-well plate. If the cells are resuspended in 5 mL of Culture Medium, there will be enough cell suspension for about 50 wells. To quickly observe osteoclasts formation, culture the cells at a higher density.
6. Feed the cells with 100  $\mu$ L of Culture Medium every 3 - 4 days. Cells will begin to fuse and form osteoclasts around day 5 (fig 1).
7. Count the osteoclasts by staining with tartrate-resistant acid phosphatase (TRAP Staining Kit, Cat. No. KT-008).





Figure 1: Osteoclasts differentiation


## EXAMPLES

1. TRAP Staining Kit (Cat. No. KT-008):  
Osteoclasts were fixed then stained with 5 mL of a mixture containing chromogenic substrate and tartrate-containing buffer.



TRAP Staining


2. Quantitation of TRAP in culture supernatant (Cat. No. KT-008):  
Thirty microliters of culture supernatant was incubated for 3 hours in the presence of chromogenic substrate/tartrate-containing buffer. The samples were read at wavelength 540 nm.



Measurement of TRAP in Osteoclasts culture supernatant

### 3. Pit Assay:

Primary precursor osteoclasts cultured on ivory for 7 – 14 days. The section was sonicated in 5 mL of 1M ammonia solution to disrupt the cells. The ivory section was stained with Mayer's hematoxylin solution for 1 minute then washed and dried.



Resorption pits on ivory section (HE staining)

### 4. Scanning electron microscopy (SEM):

SEM of the ivory section used in the Pit assay.



Resorption pits on ivory section

## **FOR RESEARCH USE ONLY**

**KAMIYA BIOMEDICAL COMPANY**

12779 Gateway Drive, Seattle WA 98168

Tel: (206) 575-8068 Fax: (206) 575-8094

Email: LifeScience@k-assay.com

www.k-assay.com