

KAMIYA BIOMEDICAL COMPANY

S-100 β ELISA

**For the quantitative determination
of S-100 β in human, mouse and rat plasma.**

Cat. No. KT-435

For Research Use Only.

PRODUCT INFORMATION

S-100 β ELISA Cat. No. KT-435

INTENDED USE

The S-100 β ELISA is for the quantitative determination of S-100 β in human, mouse and rat plasma. For research use only.

INTRODUCTION

S-100 protein has a molecular weight of 21K Dalton and consists of two subunits, α chain and β chain. It is known that combination of these subunits is different from the location in human body. S-100 $\beta\beta$ is localized in glial cell and schwann cell, S-100 $\alpha\beta$ in glial cell and S-100 $\alpha\alpha$ in striated muscle, heart and kidney.

It was reported that the concentration of S-100 β in cerebrospinal fluid was a useful marker for studying the degree of brain damage after head injury, cerebral hemorrhage and ischemic stroke. Recently another report described that the increasing of S-100 β in blood correlated to the degree of brain damage after cerebral ischemia, infarction, hemorrhage and severe head injury.

PRINCIPLE

This ELISA kit for determination of S-100 β in plasma samples is based on the sandwich enzyme immunoassay. During the first immune reaction, the S-100 β in calibrators or samples binds to the rabbit anti-bovine S-100 β antibodies which are coated on the surface of the microtiter plate. After incubation and plate washing, labeled antibody (biotinylated rabbit anti-bovine S-100 β antibodies) are added to bind to the antigen-antibody complex. Then HRP-labeled streptoavidin is added to form biotinylated rabbit anti-bovine S-100 β -antigen-antibody complexes. Finally, HRP enzyme activity is determined by o-phenylenediamine dihydrochloride (OPD) and the concentration of S-100 β is calculated.

COMPONENTS

Component	Form	Quantity	Main Ingredient
1. Antibody-Coated Plate	MTP ^{*1}	1 plate (96-well)	Rabbit anti-bovine S-100 β
2. S-100 β Calibrator	Lyophilized	1 vial (5 ng)	Bovine S-100 β
3. Labeled Antibody	Liquid	1 bottle (11 mL)	Biotinylated rabbit anti-bovine S-100 β
4. SA-HRP solution	Liquid	1 bottle (11 mL)	HRP labeled streptoavidin
5. Substrate Buffer	Liquid	1 bottle (26 mL)	Citrate buffer containing 0.015% hydrogen peroxide
6. OPD Tablet	Tablet	2 tablets	o-Phenylenediamine dihydrochloride
7. Stop Solution	Liquid	1 bottle (11 mL)	1M H ₂ SO ₄
8. Buffer Solution	Liquid	1 bottle (20 mL)	Phosphate buffer
9. Wash Solution (Concentrate)	Liquid	1 bottle (50 mL)	Concentrated saline
10. Plate Seal		4 sheets	

MTP^{*1}..... Microtiter plate

MATERIALS REQUIRED BUT NOT PROVIDED

- Photometer for microtiter plate (plate reader), which can read absorbance up to 2.5 at 490 nm
- Microtiter plate shaker
- Washing device for microtiter plate, dispenser with aspiration system
- Micropipettes, multi-channel pipettes for 8 wells or 12 wells and their tips
- Glass test tubes for preparation of Calibrator Solution
- Graduated cylinder (1,000 mL)

- Distilled water or de-ionized water

NOTES

1. EDTA-2Na (1 mg/mL) additive blood collection tubes are recommended for the plasma sample collection. It is strongly recommended that plasma samples should be used as soon as possible after collection. If the samples are tested later, they should be aliquoted and frozen at or below -30°C and thawed before the assay is run. Avoid multiple freeze / thaw cycles.
2. S-100 β calibrators and substrate solution should be prepared just before use. This kit can be used dividedly in strips of the plate. In such case, the rest of reconstituted calibrator in glass vials or tubes should be stored at 4°C (stable for 2 weeks). It is also possible to keep calibrator stable for 4 weeks if calibrator be stored at or below -30°C.
3. During storage of the concentrated wash solution at 4°C, precipitates may be observed, however they will dissolve when diluted.
4. When concentration of S-100 β in samples is expected to exceed 5 ng/mL, the sample needs to be diluted with buffer solution to a proper concentration.
5. Read optical absorbance of reaction solution in the wells immediately after stopping the color reaction.
6. Perform all the determinations in duplicate.
7. Protect reagents from strong light (e.g. direct sunlight) during storage and assay.
8. Satisfactory performance of the test is guaranteed only when reagents are used from a kit with identical lot number.
9. As pipetting operations may affect the precision of the assay, precisely pipette the prepared Calibrator Solutions or samples into corresponding wells. Use a new tip for each sample and calibrator to avoid cross-contamination. Use clean test tubes or vessels.
10. Always run a calibration curve when testing samples.
11. The total pipetting time of calibrator solutions and samples for a whole plate should not exceed 30 minutes.
12. During incubation except the color reaction, the plate should be shaken gently with a microtiter plate shaker to promote immunoreaction (approximately 100 rpm).

REAGENT PREPARATION

1. Preparation of Calibrator:
Reconstitute the calibrator (lyophilized S-100 β 5 ng/vial) with 1 mL of buffer solution, which makes a 5 ng/mL calibrator solution. The reconstituted calibrator solution is to be diluted with the same volume of buffer solution (e.g. 0.2 mL calibrator + 0.2 mL buffer solution), to make a 2.5 ng/mL calibrator solution. Repeat the dilution to make each calibrator of 1.25, 0.625, 0.313, 0.156 and 0.078 ng/mL. Buffer solution is to be used as the 0 ng/mL calibrator.

Note: Calibrator Solution must be prepared immediately before assay. Use clean test tubes or vessels.

2. Preparation of Substrate Solution: Dissolve one OPD Tablet with 12 mL of Substrate Buffer.

Note: Substrate Solution must be prepared immediately before assay. Use clean test tubes or vessels.

3. Preparation of Wash Solution: Dilute 50 mL of Wash Solution Concentrate to 1,000 mL with distilled or de-ionized water.
4. Other reagents are ready for use.

STORAGE

Store kit at 4°C.


ASSAY PROTOCOL

1. Bring all reagents and samples to room temperature (20-30°C) before starting the assay.
2. Fill 0.30 mL/well of washing solution into the wells and aspirate the solution in the wells. Repeat this washing procedure further twice (total 3 times). Finally, invert the plate and tap it onto an absorbent surface, such as paper towels, to ensure removal of most of the residual wash solution.

3. Fill 100 μL of buffer solution into all wells first and then introduce 20 μL of calibrator solutions (0, 0.078, 0.156, 0.313, 0.625, 1.25, 2.5, 5 ng/mL) or samples into the wells. The total pipetting time of calibrator solutions and samples for a whole plate should not exceed 30 minutes.
4. Cover the plate with a Plate Seal and incubate it at room temperature for 3 hours. During the incubation, the plate should be shaken with a microtiter plate shaker. (approximately 100 rpm)
5. Take off the Plate Seal, aspirate and wash the wells four times with approximately 0.3 mL/well of wash solution. Finally, invert the plate and tap it onto an absorbent surface, such as paper towels, to ensure removal of most of the residual wash solution.
6. Pipette 100 μL of Labeled Antibody into the wells.
7. Cover the plate with a Plate Seal and incubate it for one hour at room temperature. During the incubation, the plate should be shaken with a microtiter plate shaker. (approximately 100 rpm)
8. Take off the Plate Seal, aspirate the solution in the wells and wash the wells four times with approximately 0.3 mL/well of washing solution. Finally, invert the plate and tap it onto an absorbent surface, such as paper towels, to ensure removal of most of the residual wash solution.
9. Pipette 100 μL of SA-HRP solution into the wells.
10. Cover the plate with a Plate Seal and incubate it for one hour at room temperature. During the incubation, the plate should be shaken with a microtiter plate shaker. (approximately 100 rpm)
11. Dissolve one OPD tablet with 12 mL of substrate buffer. Make sure to prepare immediately before use.
12. Take off the Plate Seal, aspirate the solution in the wells and wash the wells five times with approximately 0.3 mL/well of washing solution. Finally, invert the plate and tap it onto an absorbent surface, such as paper towels, to ensure removal of most of the residual wash solution.
13. Pipette 100 μL of the substrate solution containing OPD into the wells, cover with a Plate Seal and incubate the plate for 30 minutes at room temperature for the color reaction.
14. Add 100 μL of stop solution into the wells to stop color reaction.
15. Read the optical absorbance of the wells at 490 nm. The dose-response curve of this assay fits best to a 4 (or 5)-parameter logistic equation. The results of unknown samples can be calculated with any computer program having a 4 (or 5)-parameter logistic function. Otherwise calculate mean absorbance values of wells containing calibrators and plot a calibration curve on semilogarithmic graph paper (abscissa: concentration of calibrator; ordinate: absorbance values). Use the average absorbance of each sample to determine the corresponding value by simple interpolation from this calibration curve.

PERFORMANCE

Typical Calibration Curve (example only, a new calibration curve for each run must be established by the end-user)



Analytical Recovery

<Human plasma A>

Added Bovine S-100□□ (ng/mL)	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
0.00	0.21		
0.10	0.30	0.31	96.77
0.50	0.67	0.71	94.37
2.00	1.91	2.21	86.43

<Human plasma B>

Added Bovine S-100□	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
0.00	0.18		
0.10	0.27	0.28	96.43
0.50	0.61	0.68	89.71
2.00	1.65	2.18	75.69

<Rat plasma A>

Added Bovine S-100□□	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
0.00	0.55		
0.10	0.62	0.65	95.39
0.50	0.89	1.05	84.76
2.00	1.99	2.55	78.04

<Rat plasma B>

Added Bovine S-100□	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
0.00	1.03		
0.10	1.11	1.13	98.23
0.50	1.38	1.53	90.20
2.00	2.27	3.03	74.92

<Rat plasma C>

Added Bovine S-100□	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
0.00	0.45		
0.10	0.53	0.55	96.36
0.50	0.79	0.95	83.16
2.00	1.78	2.45	72.65

<Mouse plasma A>

Added Bovine S-100□□ (ng/mL)	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
0.00	0.08		
0.10	0.16	0.18	88.89
0.50	0.53	0.58	91.38
2.00	1.81	2.08	87.02

<Mouse plasma B>

Added Bovine S-100□□ (ng/mL)	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
0.00	0.11		
0.10	0.21	0.21	100.00
0.50	0.54	0.61	88.53
2.00	1.72	2.11	81.52

Precision and reproducibility

- Intra-assay CV (%): 2.99 – 4.82
- Inter-assay CV (%): 4.82 – 9.20

Assay Range

0.078 – 5 ng/mL

Cross-Reactivity

The ELISA kit shows 0.2% cross-reactivity to bovine S-100αα.

FOR RESEARCH USE ONLY**KAMIYA BIOMEDICAL COMPANY**

12779 Gateway Drive, Seattle, WA 98168

Tel: (206) 575-8068 Fax: (206) 575-8094

Email: LifeScience@k-assay.com

www.k-assay.com