



KAMIYA BIOMEDICAL COMPANY

Mouse Alpha-1 Acid Glycoprotein ELISA

**For the quantitative determination of alpha-1 acid glycoprotein
in mouse serum or plasma.**

Cat. No. KT-685

For Research Use Only.

PRODUCT INFORMATION

Mouse Alpha-1 Acid Glycoprotein ELISA Cat. No. KT-685

PRODUCT

The **K-ASSAY®** Mouse Alpha-1 Acid Glycoprotein ELISA is an enzyme immunoassay for the quantitative determination of alpha-1 acid glycoprotein (α -1-AGP) in mouse serum or plasma. For research use only.

INTRODUCTION

α -1-AGP is an acute phase protein that is elevated in serum as a result of injury, infection or disease. Studies have demonstrated a tenfold increase in α -1-AGP in BALB/c mice after injection with lipopolysaccharide.

PRINCIPLE

The **K-ASSAY®** Mouse α -1-AGP ELISA is based on a solid phase enzyme-linked immunosorbent assay (ELISA). The assay uses affinity purified anti-mouse α -1-AGP antibodies for solid phase (microtiter wells) immobilization and horseradish peroxidase (HRP) conjugated anti-mouse α -1-AGP antibodies for detection. The test sample is diluted and incubated in the microtiter wells for 45 minutes. The microtiter wells are subsequently washed and HRP conjugate is added and incubated for 45 minutes. This results in α -1-AGP molecules being sandwiched between the immobilization and detection antibodies. The wells are then washed to remove unbound HRP-labeled antibodies and TMB Reagent is added and incubated for 20 minutes at room temperature. This results in the development of a blue color. Color development is stopped by the addition of Stop Solution, changing the color to yellow, and optical density is measured spectrophotometrically at 450 nm. The concentration of α -1-AGP is proportional to the optical density of the test sample.

COMPONENTS

- Anti-mouse α -1-AGP antibody coated microtiter plate with 96 wells (provided as 12 detachable strips of 8)
- Enzyme Conjugate Reagent, 11 mL
- Mouse α -1-AGP Calibrator (lyophilized)
- Diluent (10X), 25 mL
- Wash Solution (20X), 50 mL
- TMB Reagent (One-Step), 11 mL
- Stop Solution (1N HCl), 11 mL

MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes and tips
- Distilled or de-ionized water
- Polypropylene or glass tubes
- Vortex mixer
- Absorbent paper or paper towels
- Micro-Plate incubator/shaker mixing speed of ~150 rpm
- Plate washer
- A microtiter plate reader capable of measuring absorbance at 450 nm, with a bandwidth of 10 nm or less and an OD range of 0-4 OD
- Graph paper (PC graphing software is optional)

GENERAL INSTRUCTIONS

All reagents should be allowed to reach room temperature (18-25°C) before use.

Serum or plasma samples should be diluted ~15,000 fold with 1X diluent in order to obtain values within the calibration range.

WASH SOLUTION PREPARATION

The wash solution is provided as a 20X stock. Prior to use dilute the contents of the bottle (50 mL) with 950 mL of distilled or de-ionized water.

DILUENT PREPARATION

The diluent is provided as a 10X stock. Prior to use estimate the final volume of diluent required for your assay and dilute one (1) volume of the 10X stock with nine (9) volumes of distilled or de-ionized water.

CALIBRATOR PREPARATION

The mouse α -1-AGP calibrator is provided as a lyophilized stock.

1. Reconstitute the lyophilized mouse α -1-AGP reference calibrator by adding the volume of de-ionized or distilled water indicated on the vial label and mix gently until dissolved.
2. Label 8 polypropylene or glass tubes as 300, 150, 75, 37.5, 18.8, 9.4, 4.7 and 0 ng/mL.
3. Into the tube labeled 300 ng/mL, prepare the 300 ng/mL working calibrator as detailed on the stock vial label.
4. Dispense 250 μ L of 1X diluent into the remaining tubes.
5. Prepare the 150 ng/mL calibrator by diluting and mixing 250 μ L of the 300 ng/mL calibrator with 250 μ L of diluent in the tube labeled 150 ng/mL.
6. Similarly prepare the 75, 37.5, 18.8, 9.4 and 4.7 ng/mL calibrators by serial dilution.

Please Note: The reconstituted calibrator remains stable for at least 10 days at 4°C but should be aliquoted and stored frozen at -20°C if future use is intended.

SAMPLE PREPARATION

General Note: Our studies find that α -1-AGP is present in BALB/c mouse serum at concentrations of 0.16 to 1.4 mg/mL. In order to obtain values within the range of the calibration curve we suggest that samples initially be diluted 15,000 fold using the following procedure for each sample to be tested:

1. Dispense 998 μ L and 290 μ L of 1X diluent into separate tubes.
2. Pipette and mix 2 μ L of the serum/plasma sample into the tube containing 998 μ L of diluent. This provides a 500 fold diluted sample.
3. Mix 10 μ L of the 500 fold diluted sample with the 290 μ L of diluent in the second tube. This provides a 15,000 fold dilution of the sample.
4. Repeat this procedure for each sample to be tested.

Levels of α -1-AGP may vary with mouse strain, animal husbandry conditions and study protocol. Therefore, please be aware that optimal serum or plasma dilutions should be determined empirically.

ASSAY PROCEDURE

1. Secure the desired number of coated wells in the holder.
2. Dispense 100 μ L of calibrators and samples into the wells (we recommend that samples be tested in duplicate).
3. Incubate on an orbital micro-plate shaker at 100-150 rpm at room temperature (18-25°C) for 45 minutes.
4. Wash and empty the microtiter wells 5 times with 400 μ L 1X wash solution using a plate washer.
5. Strike the wells sharply onto adsorbent paper or paper towels to remove all residual droplets.
6. Add 100 μ L of enzyme conjugate reagent into each well.
7. Incubate on an orbital micro-plate shaker at 100-150 rpm at room temperature (18-25°C) for 45 minutes.
8. Wash as detailed in 4 and 5 above.
9. Dispense 100 μ L of TMB Reagent into each well.
10. Gently mix on an orbital micro-plate shaker at 100-150 rpm at room temperature (18-25°C) for 20 minutes.
11. Stop the reaction by adding 100 μ L of Stop Solution to each well.
12. Gently mix. *It is important to make sure that all the blue color changes to yellow.*
13. Read the optical density at 450 nm with a microtiter plate reader within 15 minutes.

CALCULATION OF RESULTS

1. Calculate the average absorbance values (A_{450}) for each set of reference calibrators, and samples.
2. Construct a calibration curve by plotting the mean absorbance obtained from each reference calibrator against its concentration in ng/mL on linear graph paper, with absorbance values on the vertical or Y-axis and concentrations on the horizontal or X-axis.
3. Using the mean absorbance value for each sample, determine the corresponding concentration of α -1-AGP in ng/mL from the calibration curve.

4. Multiply the derived concentrations by the dilution factor to determine the actual concentration of α -1-AGP in the serum/plasma sample.
5. If available, PC graphing software may be used for the above steps.
6. If the OD₄₅₀ values of samples fall outside the calibration curve when tested at the suggested dilution of 15,000, samples should be diluted appropriately and re-tested.

TYPICAL CALIBRATION CURVE

A typical calibration curve with optical density readings at 450 nm on the Y axis against α -1-AGP concentrations on the X axis is shown below. This curve is for the purpose of illustration only and should not be used to calculate unknowns. Each user should obtain his or her data and calibration curve in each experiment.

α -1-AGP (ng/mL)	Absorbance (450 nm)
300	3.717
150	2.743
75	1.489
37.5	0.987
18.8	0.579
9.4	0.408
4.7	0.339
0	0.204

STORAGE

The unused kit should be stored at 4°C and the microtiter plate should be kept in a sealed bag with desiccant to minimize exposure to damp air. Test kits will remain stable until the expiration date provided that the components are stored as described above.

FOR RESEARCH USE ONLY

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168
 Tel: (206) 575-8068 Fax: (206) 575-8094
 Email: LifeScience@k-assay.com
 www.k-assay.com