

**KAMIYA BIOMEDICAL COMPANY**

# **Hydrogen Peroxide Fluorescent Detection Kit**

**For detection in fresh urine, buffers, and tissue culture media**

**Cat. No. KT-733**

**For Research Use Only. Not for Use in Diagnostic Procedures.**

**PRODUCT INFORMATION**  
**Hydrogen Peroxide Fluorescent Detection Kit**  
**Cat. No. KT-733**

**Background**

Hydrogen peroxide was first described in 1818 by Louis Jacques Thénard. Today, industrially, hydrogen peroxide is manufactured almost exclusively by the autoxidation of a 2-alkyl-9,10-dihydroxyanthracene to the corresponding 2-alkyl anthraquinone in the Riedl-Pfleiderer or anthraquinone process.

In biological systems incomplete reduction of O<sub>2</sub> during respiration produces superoxide anion (O<sub>2</sub><sup>-</sup>), which is spontaneously or enzymatically dismutated by superoxide dismutase to H<sub>2</sub>O<sub>2</sub>. Many cells produce low levels of O<sub>2</sub><sup>-</sup> and H<sub>2</sub>O<sub>2</sub> in response to a variety of extracellular stimuli, such as cytokines (TGF-β1, TNF-α, and various interleukins), peptide growth factors (PDGF; EGF, VEGF, bFGF, and insulin), the agonists of heterotrimeric G protein-coupled receptors (GPCR) such as angiotensin II, thrombin, lysophosphatidic acid, sphingosine 1-phosphate, histamine, and bradykinin, and by shear stress. The addition of exogenous H<sub>2</sub>O<sub>2</sub> or the intracellular production in response to receptor stimulation affects the function of various proteins, including protein kinases, protein phosphatases, transcription factors, phospholipases, ion channels, and G proteins. In 1894, Fenton described the oxidation of tartaric acid by Fe<sup>2+</sup> and H<sub>2</sub>O<sub>2</sub>. H<sub>2</sub>O<sub>2</sub> and O<sub>2</sub> may participate in the production of singlet oxygen and peroxynitrite and the generation of these species may be concurrent with reactions involving iron, and under some circumstances they might be important contributors to H<sub>2</sub>O<sub>2</sub> toxicity.

A substantial portion of H<sub>2</sub>O<sub>2</sub> lethality involves DNA damage by oxidants generated from iron-mediated Fenton reactions. Damage by Fenton oxidants occurs at the DNA bases or at the sugar residues. Sugar damage is initiated by hydrogen abstraction from one of the deoxyribose carbons, and the predominant consequence is eventual strand breakage and base release.

**PRINCIPLE**

The **K-ASSAY**® Hydrogen Peroxide Fluorescent Detection Kit is designed to quantitatively measure H<sub>2</sub>O<sub>2</sub> in a variety of samples. Please read the complete kit insert before performing this assay. A hydrogen peroxide calibrator is provided to generate a calibration curve for the assay and all samples should be read off the calibration curve. Samples are mixed with the Fluorescent Substrate and the reaction initiated by addition of horseradish peroxidase. The reaction is incubated at room temperature for 15 minutes. The HRP reacts with the substrate in the presence of hydrogen peroxide to convert the colorless substrate into a fluorescent product. The fluorescent product is read at 590 nm with excitation at 570 nm. Increasing levels of H<sub>2</sub>O<sub>2</sub> cause a linear increase in fluorescent product.

**COMPONENTS**

- |                                      |          |
|--------------------------------------|----------|
| • Black 96 Well Half Area Plates     | 2 Plates |
| • Hydrogen Peroxide Calibrator       | 220 µL   |
| • Fluorescent Detection Reagent      | 5 mL     |
| • Assay Buffer Concentrate           | 25 mL    |
| • Horseradish Peroxidase Concentrate | 60 µL    |

**Storage**

All components of this kit should be stored at 4°C until the expiration date of the kit.

**Other Materials Required**

- Repeater pipet with disposable tips capable of dispensing 25 µL.
- 96 well plate reader capable of reading fluorescence at 580-590 nm with excitation at 570-580 nm. Set plate parameters for a 96-well Corning Costar 3694 plate.
- Software for converting fluorescent intensity readings from the plate reader and carrying out four parameter logistic curve (4PLC) fitting. Contact your plate reader manufacturer for details.

## Precautions

As with all such products, this kit should only be used by qualified personnel who have had laboratory safety instruction. The complete insert should be read and understood before attempting to use the product. The supplied hydrogen peroxide calibrator contains very dilute H<sub>2</sub>O<sub>2</sub>.

## Sample Types

Samples that need to be stored after collection should be stored at -70°C or lower, preferably after being frozen in liquid nitrogen. Urine samples can be used after being diluted  $\geq$  1:10. This assay has been validated for buffer and media samples.

## Sample Preparation

Dilute samples  $\geq$  1:10 with Assay Buffer prior to running in the assay.

## Reagent Preparation

### Assay Buffer

Dilute Assay Buffer Concentrate 1:5 by adding one part of the concentrate to four parts of deionized water. Once diluted this is stable at 4°C for 3 months.

### Horseradish Peroxidase (HRP) Preparation

Dilute the HRP Stock solution 1:100 with Assay Buffer using the table below:


HRP Dilution Table

	1/2 Plate	One Plate	Two Plates
HRP Stock	15 $\mu$ L	30 $\mu$ L	55 $\mu$ L
Assay Buffer	1.485 mL	2.97 mL	5.445 mL
Total Volume	1.5 mL	3 mL	5.5 mL

### Calibrator Preparation

Hydrogen Peroxide Calibrators are prepared by labeling seven tubes as #1 through #7. Briefly vortex to mix the vial of H<sub>2</sub>O<sub>2</sub> calibrator. Pipet 450  $\mu$ L of Assay Buffer into tube #1 and 200  $\mu$ L into tubes #2 to #7. Carefully add 50  $\mu$ L of the H<sub>2</sub>O<sub>2</sub> Calibrator to tube #1 and vortex completely. Take 200  $\mu$ L of the solution in tube #1 and add it to tube #2 and vortex completely. Repeat this for tubes #3 through #7. The concentration of H<sub>2</sub>O<sub>2</sub> in tubes 1 through 7 will be 10, 5, 2.5, 1.25, 0.625, 0.313 and 0.156  $\mu$ M.

Use all Calibrators within 2 hours of preparation.



	Std 1	Std 2	Std 3	Std 4	Std 5	Std 6	Std 7
Assay Buffer ( $\mu$ L)	450	200	200	200	200	200	200
Addition	Stock	Std 1	Std 2	Std 3	Std 4	Std 5	Std 6
Vol of Addition ( $\mu$ L)	50	200	200	200	200	200	200
Final Conc ( $\mu$ M)	10	5	2.5	1.25	0.625	0.313	0.1569

## ASSAY PROTOCOL

1. Use the plate layout sheet on the last page of the insert to aid in proper sample and calibrator identification. Set plate parameters for a 96-well Corning Costar 3694 plate.
2. Pipet 50  $\mu$ L of samples or appropriate calibrators into duplicate wells in the plate.

3. Pipet 50  $\mu\text{L}$  of Assay Buffer into duplicate wells as the Zero calibrator.
4. Add 25  $\mu\text{L}$  of Fluorescent Substrate to each well using a repeater pipet.
5. Initiate the reaction by adding 25  $\mu\text{L}$  of the HRP Preparation to each well using a repeater pipet.
6. Incubate at room temperature for 15 minutes.
7. Read the fluorescent emission at  $585 \pm 5$  nm with excitation at  $575 \pm 5$  nm. Please contact your plate reader manufacturer for suitable filter sets.

## Calculation of Results

Average the duplicate FLU readings for each calibrator and sample. Create a calibration curve by reducing the data using computer software capable of generating a four-parameter logistic curve (4PLC) fit, after subtracting the mean FLUs for the Zero wells. The sample concentrations obtained should be multiplied by the dilution factor to obtain neat sample values.

## Typical Data


Always run your own calibration curve for calculation of results.

Sample	Mean FLU	Net FLU	H <sub>2</sub> O <sub>2</sub> Conc. ( $\mu\text{M}$ )
Zero	3,782	0	0
Standard 1	36,417	32,635	10
Standard 2	21,919	18,137	5
Standard 3	13,134	9,352	2.5
Standard 4	8,333	4,551	1.25
Standard 5	6,072	2,290	0.625
Standard 6	5,031	1,249	0.313
Standard 7	4,398	616	0.156
Sample 1	6,578	2,796	0.76
Sample 2	24,680	20,898	5.85

Do not use this data.

**Conversion Factor: 100 nM of Hydrogen Peroxide is equivalent to 3.4 ng/mL.**

## Typical Calibration Curve



## Sensitivity and Limit of Detection

Sensitivity was calculated by comparing the FLUs for twenty wells run for each of the zero and calibrator #7. The detection limit was determined at two (2) standard deviations from the zero along the calibration curve.

**Sensitivity was determined as 0.038 µM. This is equivalent to 1.9 pmol (64.6 pg) H<sub>2</sub>O<sub>2</sub> per well**

The Limit of Detection was determined in a similar manner by comparing the FLUs for twenty wells run for each of the zero and a low concentration human sample.

**The Limit of Detection was determined as 0.052 µM. This is equivalent to 2.6 pmol (88.4 pg) H<sub>2</sub>O<sub>2</sub> per well**

## Intra Assay Precision

Three buffer samples were run in replicates of 20 in an assay. The mean and precision of the calculated concentrations were:

Sample	H <sub>2</sub> O <sub>2</sub> Conc. (µM)	%CV
1	6.27	3.6
2	3.21	3.8
3	0.98	5.7

## Inter Assay Precision


Three buffer samples were run in duplicates in fourteen assays run over multiple days by three operators. The mean and precision of the calculated concentrations were:

Sample	H <sub>2</sub> O <sub>2</sub> Conc. (µM)	%CV
1	5.86	4.3
2	3.00	7.0
3	0.88	12.1

## Linearity

Linearity was determined by taking two RPMI-1640 media samples with known H<sub>2</sub>O<sub>2</sub> concentrations and mixing them in the ratios given below. The measured concentrations were compared to the expected values based on the ratios used.

High RPMI	Low RPMI	Observed Conc. ( $\mu\text{M}$ )	Expected Conc. ( $\mu\text{M}$ )	% Recovery
80%	20%	1.23	1.32	93.5
60%	40%	1.12	1.20	93.2
40%	60%	0.97	1.09	89.2
20%	80%	0.93	0.97	95.5
			<b>Mean Recovery</b>	<b>92.8%</b>



	A	B	C	D	E	F	G	H
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								

**FOR RESEARCH USE ONLY**

**KAMIYA BIOMEDICAL COMPANY**

12779 Gateway Drive, Seattle WA 98168  
Tel: (206) 575-8068 Fax: (206) 575-8094  
Email: LifeScience@k-assay.com  
www.k-assay.com