



PRODUCT INFORMATION

Copper Assay (3,5-DiBr-PAESA Chromogenic Method)

Cat. No. KT-758 / KT-759

INTENDED USE

The Copper Assay (3,5-DiBr-PAESA Chromogenic Method) is a direct colorimetric assay for the quantitative determination of copper in biological samples that does not require sample deproteinization. For research use only. Not for use in diagnostic procedures in the U.S.

PRINCIPLE

Copper from ceruloplasmin-copper complex is dissociated by a weak acidic buffer and reduced by ascorbic acid. The Cu^+ ions form a blue chelate with the 3,5-DiBr-PAESA chromogen. The intensity of this colored complex is proportional to the copper concentration in the sample and is measured at 580 nm.

COMPONENTS

	100 tests	200 tests
	<u>KT-758</u>	<u>KT-759</u>
1. Buffer, 24 mL	x1	x2
2. Chromogen (3,5-DiBr-PAESA), 0.5 mL	x1	x2
4. Copper Calibrator (200 $\mu\text{g}/\text{dL}$), 1.2 mL	x1	x2

Store all kit components at 4°C.

PRECAUTIONS

1. Fluctuating incubation temperature may result in variable results.
2. Use disposable test tubes and glassware washed with 1M HNO_3 or 1M HCl solution and distilled water.

3. Sample and reagent pipetting accuracy may affect assay performance. Please note that samples, calibrator, and Working Solution must be dispensed accurately at the μL level.
4. The temperature of the reaction may affect the O.D. reading. Please extend or shorten the chromogen reaction time depending on the ambient room temperature if necessary.
5. For cell lysates or the tissue extraction samples, a high concentration of proteins or lipids may affect the assay result. For best results, remove proteins or lipids by ultrafiltration or centrifugation.
6. Heme-containing copper cannot be measured by this assay kit.

SAMPLE PREPARATION

1. Serum or plasma

Insoluble substances in serum and plasma samples should be removed by filtration or centrifugation. EDTA-plasma samples cannot be used as EDTA interferes with this assay.

2. Tissue extracts, cell lysates, and other samples such as urine or other biological fluids:

If the sample is turbid, centrifuge at 6,000 rpm for 15 min. Collect the supernatant and use for the assay

If necessary, add small amounts of 6M HCl to the sample and adjust pH to 2.0 - 3.0. For example, add ~5-10 μL of 6M HCl per 1 mL of sample.

3. Tissue samples

Add 5% TCA solution, vortex 1 min. and incubate at 4 - 8°C for 30 min. Centrifuge at 6,000 rpm for 15 min. Collect the supernatant and use for the assay.

Sample pH should be between pH 2.0 to pH 8.0.

REAGENT PREPARATION

1. Prepare enough Working Solution for your experiment:

Working Solution			
	1 test	50 tests	100 tests
Buffer	240 µL	12 mL	24 mL
Chromogen	5 µL	250 µL	500 µL

Working Solution should be stored at 4°C and used within one month after preparation.

2. Copper Calibrator is ready to use.
3. Bring all reagents to room temperature before use.

ASSAY PROTOCOL (Microplate and Microplate Reader)

(Total reaction volume = 252 µL)

1. Add 12 µL of Blank (purified water), Copper Calibrator, or Sample to each well.
2. Add 240 µL of Working Solution to each well, mix, and incubate at room temperature for 10 minutes. Mix carefully using a pipette to avoid foaming. If a plate mixer is used for mixing, there is a risk of obtaining poor reproducibility.
3. Read the OD absorbance at 580 nm (main) and 750 nm (sub). Acceptable wavelength range: 570 - 590 nm (main) and 700 - 800 nm (sub).

Assay Protocol				
Step	(µL)	Blank	Calibrator	Sample
1	Purified water	12	-	-
	Copper Calibrator	-	12	-
	Sample	-	-	12
2	Working Solution	240	240	240
	Mix and incubate for 10 minutes at room temp.			
3	Read the OD absorbance at 580 nm (main) and 750 nm (sub).			

CALCULATION OF SAMPLE CONCENTRATION

$$\frac{(\text{OD}_{580} \text{ sample} - \text{OD}_{750} \text{ sample}) - (\text{OD}_{580} \text{ blank} - \text{OD}_{750} \text{ blank})}{(\text{OD}_{580} \text{ calib.} - \text{OD}_{750} \text{ calib.}) - (\text{OD}_{580} \text{ blank} - \text{OD}_{750} \text{ blank})} \times 200 = \text{Copper } (\mu\text{g/dL})$$

Unit Conversion:

$$\text{Copper } (\mu\text{g/dL}) \times 0.1575 = \text{Copper } (\mu\text{M})$$

Assay Example

	OD (580 nm)	OD (750 nm)	OD	ΔOD	Copper (µg/dL)
Blank	0.069	0.028	0.041	-	-
Calibrator	0.160	0.054	0.106	0.065	-
Sample	0.106	0.038	0.068	0.027	83.1

When assaying diluted samples, multiply the result by the dilution factor.

PERFORMANCE

Assay Range: 3 - 400 µg/dL

Precision: Precision was evaluated using commercially available quality control serum.

Within Run Precision	Mean (µg/dL)	S.D.	C.V.%
Level 1	76.8	2.2	2.9
Level 2	178.6	5.4	3.0

Interference:

Conjugated bilirubin	No interference up to at least 40 mg/dL
Unconjugated bilirubin	No interference up to at least 40 mg/dL
Hemoglobin	No interference up to at least 0.1 g/dL
Chyle	No interference up to at least 500 FTU

Shelf life: Until expiration date at 4°C. After opening any of the kit components, store at 4°C and use within one month. **Do not freeze.**

FOR RESEARCH USE ONLY.

Not for use in diagnostic procedures in the U.S.

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168

Tel: (206) 575-8068 Fax: (206) 575-8094

Email: LifeScience@k-assay.com

www.k-assay.com